

PRODUCT INFORMATION SHEET

DESCRIPTION:

Q-SOL is a Chlorinated Disinfectant/ Sanitiser that can be used as a non-rinse sanitiser in food preparation areas or as a hospital grade disinfectant for outbreak cleaning or blood spills. Its convenient sachet packaging allows the right concentration of Chlorine to be used every time.

FEATURES AND BENEFITS:

- Chlorinated hospital grade disinfectant
- Convenient sachet packaging
- Easy use as non-rinse sanitiser
- Consistent chlorine content

TECHNICAL INFORMATION:

Appearance	Pink Powder
Odour	Slight Chlorine
Specific Gravity	1.00
Solubility in Water	20g per 100g of water at 20 °C

USE:

Thoroughly pre-clean area to be sanitized or disinfected, allow minimum 10 minutes contact time if rinsing.

Non-rinse Sanitising (100 PPM):	Add 1 sachet to 20L water
Surface Sanitising (500 PPM):	Add 1 sachet to 4L water
Blood Spills (10,000 PPM):	Add 1 sachet to 200ml water
Infectious Cleaning (1,000 PPM):	Add 1 sachet to 2L water

Note: The pink grade of Q-SOL acts as a useful indicator as colour changes from pink to grey when bacteria persist in the system or equipment. If the solution remains pink at the end of the circulation or cleaning process this means sanitization is complete.

SAFETY: Contact with acids liberates toxic gas. Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation.

Keep out of reach of children. Read label before use. Avoid breathing dust, fume, gas, mist, vapours or spray. Wash hands, face and all exposed skin thoroughly after handling. Use only outdoors or in a well-ventilated area. Wear protective clothing, gloves, eye / face protection and suitable respirator.